

ESPECIFICACIONES BK2001-3001

Product Name:	Aneroid Sphygmomanometer with Stethoscope	Bladder:	PVC/Latex bladder with two tube,size 22X12cm
Model Number:	BK2001-3001	Carrying Bag:	069 leather bag
Measure scope:	0-300mmHg	Stethoscope:	Single head stethoscope
Accuracy:	+/-3mmHg	Chest piece:	Aluminum chest piece
Sub-division:	2mmHg	Y-tube:	60cm long Y-tube
Aneroid Gauge:	110g/120g/140g	Color:	black,blue,grey,red,green etc
Cuff:	023 Nylon/cotton cuff with/without D-ring 48X14cm/51X14cm		