

Overall Length	1075 mm / 42.3"
Overall Width	650 mm / 26" (Armrest Width)
Overall Height	900 mm / 35.4"
Front Wheels	190 mm / 7.4"
Middle Wheels	N/A
Rear Wheels	320 mm / 13"
Weight W / Batteries	70.9 kg / 156.3 lbs
Max. Speed	8 kmph / 5 mph
Weight Capacity	127 kg / 280 lbs
Ground Clearance	46 mm / 1.8"
Grade Climbable	10 degree
Curb Climbable	43 mm / 1.7"
Turning Radius	980 mm / 38.6"
Suspension	N/A
Brake	Electro-Mechanical
Seat Type	Padded Breathable Ballistic Nylon Back Rest W/ Adjustable Back Angle & Seat Belt
Seat Depth	440 mm / 17.3"
Armrest Height	230 mm / 9"
Seat Back Angle	93 degree
Seat Back Height	440 mm / 17.3"
Seat Width	457 mm / 18"
Drive Train	2 Motor, Rear-Wheel Drive
Motor Size	420W, 4600 r.p.m
Battery Size	(2) 12V. 36Ah
Weight of Battery	25.5 kg / 56.2 lbs
Travel Range	19 km / 12 Miles
Battery Charger	5A Off Board