### GOT YOU

## BRAND IS OUR ETERNAL PURSUIT SKB041-3

#### MEDICAL CARE SPECIALIST



## Our pricing, no middlemen make a difference











Factory Network

Customer



### $0-75^{\circ}$

#### Back-rest lifting

Allows the resident to get out of bed with confidence



 $0-30^{\circ}$ 

#### Leg-rest lifting

Protection and freedom of movement







### $600 - 920_{mm}$

#### Height Adjustment

All four side rails are raised entirely to provide patient security



## $0-18^{\circ}$

#### Trendelengurg

Protection and freedom of movement





#### Detail display



Fifth Wheel



Height adjustment panel Light and convenient



Manual operate the back rest position



Mattress
Convenient transfer of patients



Central locking panel Prevention of trolley movement



Manual operate the knee rest position



Trendelenburg panel Light and convenient



IV pole with 4 hooks Convenient to hang water



8 inch wheels Deluxe mute wheel

#### The shape is similar and the quality is different

The peers have been imitating all the time, but we have never been surpassed.

Similar in shape and different in process

### 1 Line spray exquisite workmanship

High hardness, high wear resistance






Manual spray, rough technology, easy to fall off the surface.

Similar in shape and different in process

Sturdy, durable, soft treatment

2 Fold able aluminum alloy guardrail





Cheap accessories, easy to break, and have a safety hazard.

#### Introduction

Dimension (L\*W\*H): 2100\*750\*600-920mm

Reversed trendelenburg: 0-18°

Leg—rest lifting: 0—30°

Loading bearing: 220kg

Equipped with two USA hydraulic pumps

X-ray translucent platform

Back rest lifting:  $0-75^{\circ}$ 

Trendelenburg: 0-18°

Height adjustment: 600-920mm

Electrostatic spray steel frame and ABS base

Fold able aluminum alloy guardrail

5th wheel to control the direction of moving

Pedal control elevation trendelenburg/reverse trendelenburg, Located on both sides of the table

Oxygen cylinder holder and tray for placement of x—rays cassette

I piece SS telescoping I V stand and previsions

60mm waterproof PU mattress

Rubber bumpers on the four corners

8 inch central controlled caster

# Multiple program selection, full automatic intelligent care

Reliable quality and good faith service



#### About us

Medical sickbed expert

Jiangsu Saikang Medical Equipment Co.,Ltd. has been established since 2002. It is located in Zhangjiagang City, along the Economic circle of Yangtse River Delta, enjoying the highly developed economic expansion and convenient transportation access.

#### **FDA**



#### CE

